

Sacro  Tech™

If You Thought that Your Sacroiliac Joint Belt Didn't Work, you were probably right...



The Problem:

For years doctors of Chiropractic, Physical therapists, Orthotists, and Physiatrists have been searching for a way to stabilize the Sacroiliac Joints after an injury or during what is termed an acute flair-up (an increase of inflammation and tenderness to a past trouble spot). The problem was that no S.I. Joint belts on the market could offer enough direct compression to allow the S.I. Joints to stabilize during the healing process.

The Solution:

Dr. Rick Swartzburg, D.C. utilized his clinical experience in treating Sacroiliac Joint Dysfunction and injury patients and merged it with Biocybernetics patented mechanical pulley system to bring you a brand-new belt that can actually apply enough compression to the S.I. Joints, so that they will hold their alignment and not continue to cause further injury during the healing process.

The unique Patented design will apply a customizable amount of pressure to your sacroiliac region. Simply pull on the side strap to allow the belt to cinch together, thereby applying desirable supportive compression to the sacroiliac joints. The single pull strap uses a mechanical advantage pulley system to give you several times the amount of pressure that you could manually supply to help relieve sacroiliac pain. Once you have added the correct amount of pressure for you, place the pull-tab on the belt.



One reason that doctors, therapists and their patients cannot seem to find a Sacro-Iliac Joint belt that works well is because most sacroiliac belts do not add a significant amount of pressure to the Sacroiliac joints without over compressing the whole hips and pelvis. When

you over compress the whole of the hips and pelvis, you bring an undesirable consequence of constricting the blood supply. This belt adds pressure to the S.I. joints by using a mechanical advantage system that targets the area between and directly underneath the rear pads of the belt. It does all of this with an easy, one-handed adjustment.

"Use of the SacroTech Biocybernetics Sacroiliac Joint Belt can accomplish what no other Sacroiliac Belt on the market has previously been able to do by actually offering real support to truly lock the S.I. Joints in place during its healing process. This is very necessary after a sacroiliac injury or when the Sacroiliac joints have become inflamed, as any further dysfunctional movement will cause a shearing force and further damage the joints. I have personally seen patients heal faster with use of this belt and noticed that they do not become dependent on the belt once the S.I. joints become stabilized properly over time." Swartzburg, D.C., Co-Creator of the SacroTech™